

## Welch Allyn

# Invitation

You are cordially invited to Welch Allyn's booth at

The 31th Asia Pacific Academy of Ophthalmology Congress

in conjunction with

The 57th Annual Meeting of the Ophthalmological Society of Chinese Taipei

Come and experience the latest handheld Spot™ Vision Screener

Venue: Booth A1001

Taipei World Trade Center Exhibition Hall 1

For more information, please contact:

Esther Yeo: esther.yeo@welchallyn.com

Miranda Ma: miranda.ma@welchallyn.com



#### Would you like to win one of the Welch Allyn Ophthalmoscope?

Visit our booth and complete a survey to participate in our daily lucky draw from 24<sup>th</sup> to 26<sup>th</sup> March







- For 6 month-old babies through to adults, and with limited mobility
- User friendly with touch screen display
- Successful capture 97% of time
- Easy to understand result reports
- WIFI connectivity

#### PanOptic™ Ophthalmoscope

- · Fast, easy entry into undilated pupils.
- 5x larger view of the fundus vs. standard ophthalmoscopes in an undilated eye.
- 25º field of view vs. the standard 5º field of view.





### 3.5V Coaxial Ophthalmoscope

- · Halogen HPX® lamp for true tissue color
- · Detect corneal abrasions with cobalt blue filter
- Polarizing filter virtually eliminates corneal reflection
- · Red-free filter may be used with any pressure

#### Pocket LED Ophthalmoscope

- SureColor™ LED lamp technology, provide the exact light for fundoscopic examination
- Durable pocket-style ophthalmoscope based on rigorous drop test vs other ophthalmoscopes

